



- C Added radii to corners
- B Material of Dowel Pin Spec
- A Generation Update
- 0 Original Release

8/4/22
 5/4/22
 4/28/22
 1/27/22

| | | | | | |
|---|---------------------------|---|---------|--|---------|
| UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES | | CAD GENERATED DRAWING, DO NOT MANUALLY UPDATE | | Magnitude Instruments, LLC 224 Technology Center 200 Innovation Boulevard State College, PA 16803 | |
| TOLERANCES: | | DO NOT SCALE DRAWING | | TITLE | |
| FRACTIONS | DECIMALS | ANGLES | | enclosure, enViSion, Gen II Sample Mount Plate | |
| $\pm 1/64$ | .XX $\pm .010$ | $\pm 1^\circ$ | | SIZE | SCALE |
| | .XXX $\pm .005$ | | | A | 1:2 |
| MATERIAL | aluminum .25 | DRAWN | TS | SHEET | DWG. # |
| | | | 1/27/22 | 1 | 3900132 |
| FINISH | black powder coat RAL9005 | CHECKED | TG | 1 | OF 1 |
| | | | 1/27/22 | | |

THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF MAGNITUDE INSTRUMENTS, LLC. ANY REPRODUCTION IN PART OR WHOLE WITHOUT THE WRITTEN PERMISSION OF MAGNITUDE INSTRUMENTS, LLC IS PROHIBITED.